

QUICK REFERENCE GUIDE

XLR8 · XLR8+ · UNO+

FEATURES

Therapy Mode

“C” = Continuous
“V” = Variable
“I” = Intermittent

Treatment Time

Indicates the therapy has been in progress

Target Pressure

Time and Date

Keypad Lock/Unlock Status

Red = Locked
Green = Unlocked

On/Off Button

Lock/Unlock Button

Power Jack Inlet Use for charging

Alert Notification

Icon appears whenever an event has occurred (Leak Alert, Check Dressing, Blockage, Canister Full or Low Battery).

The Events can be acknowledged by pressing Menu - Alert setup - Alert log, select each alert and press menu to clear.

Battery Indicator

Suction Strength

Green Bars = Seal is adequate
Red Bars = Seal may be compromised

System Status

Green Check (✓) = proper function
Red Exclamation (!) = problem
! = Potential small leak on dressing.

Menu/Select Button

Press to silence or lock in user changes.

Exit/Cancel Button

Press to go back to previous, or main screen.

Up/Down Arrow Buttons

Press in main screen to increase/decrease pressure setting. In Menu/Select screen, press to scroll through menu and implement functions.

QUICK REFERENCE GUIDE

XLR8 · XLR8+

Power On/Off



Press **On/Off Button**

Note: Pump will not power off if in **Locked** Mode.

Adjust Pressure Setting



Press **Up Arrow** to increase the vacuum pressure to adjust in 1 mm/hg increments.

Press **Down Arrow** to decrease the vacuum pressure to adjust in 1 mm/hg increments.

Press and Hold **Up/Down Arrow** to increase/decrease in 10 mm/hg increments.

Setting Variable Pressure



Press **Menu/Select**

Press **Up/Down** to highlight Treatment Mode then press **Menu/Select**



Press **Up/Down** to highlight Variable then press **Menu/Select**

Press **Up/Down** to select Preset Treatment times of 125/30, 80/30 and 125/0

Press **Menu/Select** to lock in

For Custom Variable Pressure Setting



SET VARIABLE HIGH PRESSURE TIME:

Press **Up/Down** in 1 minute increments. This can be set from 1 to 15 minutes. To lock in time, press **Menu/Select**



SET VARIABLE HIGH PRESSURE:

Press **Up/Down** to set high pressure. To lock in high pressure press **Menu/Select**

SET VARIABLE LOW PRESSURE TIME:

Press **Up/Down** in 1 minute increments. This can be set from 1 to 15 minutes. To lock in time, press **Menu/Select**

SET VARIABLE LOW PRESSURE:

Press **Up/Down** to set low pressure. To lock in low pressure, press **Menu/Select**

Silence Alert



Press **Menu/Select**

Note: The alert will sound again after 5 minutes if problem is not resolved.

Reset Treatment Time



Press and HOLD **Menu/Select**



Press and RELEASE **On/Off**



Press and RELEASE **Lock/Unlock**

RELEASE the **Menu/Select**

Lock/Unlock Pump



Press and hold **Lock/Unlock** and **Exit/Cancel** at the same time for 5 seconds.



When unlocked, a **green open lock** icon will appear on the bottom left of the LCD screen.

When locked, a **red closed lock** icon will appear on the bottom left of the LCD screen.

A **locked** screen disables the user from making any changes to pump settings.

Setting Intensity Mode



Press **Menu/Select**



Press **Up or Down Arrow** and highlight Intensity

Press **Menu/Select**

Press **Up or Down Arrow** to set intensity to 1, 2 or 3 (low to high)

Press **Menu/Select** to lock in the intensity setting

Clear Alert Log



User MUST be in the main screen



Press and HOLD **Menu/Select**

Press and RELEASE **Lock/Unlock**



Press and RELEASE **Exit/Cancel**

RELEASE **Menu/Select**

DMR-03-093 Rev E