

GENADYNE

CLINICAL
DECISION
GUIDE

the power to *heal*

DMR-06-097 Rev B



Disposable/Reusable	Disposable (97607 and 97608 Application Code)				Reusable (E2402 Rental Code)				Reusable (E2402 Rental Code)			
Wound Depth	1-2 cm				1-2 cm use dressing alone, 2-5 cm use cavity filler				Any depth			
Type of Wound/Dressing	All wound types • Chronic • Acute • Traumatic • Subacute and dehisced wounds • Partial-thickness burns • Ulcers (such as diabetic or pressure) • Flaps and grafts • Closed Surgical Incisions • Drainage-low				All wound types • Chronic • Acute • Traumatic • Subacute and dehisced wounds • Partial-thickness burns • Ulcers (such as diabetic or pressure) • Flaps and grafts • Closed Surgical Incisions • Drainage-low, med				All wound types • Chronic • Acute • Traumatic • Subacute and dehisced wounds • Partial-thickness burns • Ulcers (such as diabetic or pressure) • Flaps and grafts • Closed Surgical Incisions • Drainage-low, med			
Exudate	Low; dressing will absorb up to 4x it's weight in addition to the 70 ml canister				Low to medium - less than 200-300 mls x 3 days. Canisters; 2 sizes available. If exceeds more than 2 canisters weekly, upgrade to XLR8+				Low to medium, high Canisters; 5 sizes available, may be left up to one week			
Dressing Type	Peel and apply white polymer foam dressing available in 10 sizes				Hybrid peel and apply green foam dressing in 4 sizes or small and medium green foam - sm/med, PVA foam				Hybrid foam dressing; 4 sizes available, green foam thin and thick in small, medium, large and extra large, spiral, bridge, white PVA, abdominal			
Accessories	Y connector, carry bag, power cord				Y connector, silver contact layer, carry bag, power cord				Y connector, silver contact layer, carry bag, power cord, bed hanger, iv pole holder			
Dressing Change Frequency	• Peel and apply white polymer foam dressing - up to 7 days if cavity is more than 1-2cm add white PVA or green foam filler				• Hybrid peel and apply green foam dressing - up to 7 days • Green Foam dressing - 3x per week on most wounds or per physician order • If billing the Hybrid Dressing in the homecare setting under code A6550, follow Medicare guidelines.				• Hybrid peel and apply green foam dressing - up to 7 days • Green Foam dressing - 3x per week on most wounds or per physician order • If billing the Hybrid Dressing in the homecare setting under code A6550, follow Medicare guidelines.			
Dressing Description	• Polyurethane Film • Superabsorbant foam layer • Silicone gel adhesive layer				• Polyurethane Film - allows effective management of exudate away from wound into canister • Genadyne Green Foam - allows for filtering of exudate • Silicone Hybrid Dressing - allows for strong adhesion to wound site, conforms to wound shape and scar line, reduces risk of maceration to periwound				• A fully reticulated polyurethane green foam made with polyether resin to allow better hydrolytic stability, providing increased performance while under pressure • Genadyne Green Foam - allows for filtering of exudate • Silicone Hybrid Dressing - allows for strong adhesion to wound site, conforms to wound shape and scar line, reduces risk of maceration to periwound			
Tunneling/Undermining	No				No undermining, small tunneling ok with white pre-moistened PVA				Yes use white pre-moistened PVA			