

GENADYNE	UNO Plus Cleaning Instruction	DMR-07-107	
	Owner: QA/RA	Rev. A	Page 1 of 2

REVISION HISTORY			
REV	DESCRIPTION OF CHANGES	DATE	BY
1	Draft	03/19/20	AG
A	Formal Release	03/28/20	AG

GENADYNE	UNO Plus Cleaning Instruction	DMR-07-107	
	Owner: QA/RA	Rev. A	Page 2 of 2

CAUTIONARY NOTE:

- Disposable gloves, mask, and lab coat should be worn at all time.
- When disinfecting the units, extreme caution must be taken to ensure that no excess fluid is allowed to pool in the areas surrounding the screen or the outside edges of the rubber seal.
- Excessive pressure should not be applied to the screen.
- Be very careful not to allow any moisture entering the UNO+ NPWT device while cleaning.
- Do not use abrasive wipes on the screen.

Cleaning Process

1. If there is a Carrying Bag, please take it off from the device. This bag is reusable only if laundry is done.
2. Use Sporidicin® Disinfectant Towelettes/Wipes or any FDA approved disinfectant cleaning solutions/ wipes. If unsure, a Genadyne Bio-Med representative can help the facility determine the usability of the specific disinfectant. If necessary, wring out the wipe before using to reduce liquid in wipe.
3. Using the wipe, start at the canister side of the unit, wiping down the unit. Wipe down the sides of the unit and then wipe down the front of the unit. Gently clean smudge marks off the area around the screen. Tilt the unit so fluid will not pool onto screen area. See photos:
4. Wipe down the screen. Extreme caution must be taken when wiping down the screen to not allow fluid ingress into the screen. It should be done very lightly as to not damage the screen.
5. Lay unit flat with screen area up and allow to air dry.
6. After the unit has dried, wipe down the screen with “micro-fiber” cloth to reduce liquid collection on screen and clean off any smudge marks on the screen.
7. Look for obvious damage to the housing (crack/broken), screen, and power port. If there is damage then tag it and follow local procedures for non-conforming products or contact Genadyne for a RMA.
8. Use wipe to clean adapter, cord, and crevices of the power connector. Wrap the power cord in a looping fashion.
9. When the cleaning process is done, dispose the gloves and wipes. Please be sure to wash your hands.
10. Move clean unit to another area to avoid cross-contamination. Do not handle clean products with dirty gloves.