

GENADYNE

XLR8+



The power to *Heal*
in the palm of your hand.

WHEN IT COMES TO FEATURES, WE'VE THOUGHT OF EVERYTHING...

- Lightweight (1.65lbs or 750grams)
- Whisper-quiet operation
- 12+ hours of battery usage with 3 hour recharge time
- Supports multiple languages
- 5 advanced alerts - leak, blockage, canister full, low battery and critical battery
- Pump screen provides alerts to caregivers ensuring timely responses
- Simple lockable controls
- Comes with carry bag and bed hanger
- FDA approved for both Hospital and Home Care
- 3 intensity settings allowing user to choose air flow rate based on wound size and volume
- Continuous, Intermittent and Variable therapy options

Benefit of Variable Pressure

- Completely customizable high and low pressure settings
- Always maintains negative pressure at wound site
- Reduces discomfort caused by pressure spikes



GENADYNE

SUPPLIES AND ACCESSORIES

XLR8 and XLR8+

Dressing Components




Green Foam

- A fully reticulated polyurethane foam made with polyether resin to allow better hydrolytic stability, providing increased performance while under pressure
- Duo-decahedron structure gives the cells three-dimensional skeletal strands providing enhanced filtering qualities
- Promotes perfusion and assists in granulation formation
- Fraying is significantly reduced, preventing stray pieces from being left in the wound bed
- Green Foam is easier to identify in the wound bed



Port Pad

- Soft patient interface which enhances patient comfort and helps avoid possibility of pressure points
- Film-based secure system for easy application and dressing security



Transparent Film

- Easy 1-2-3 application guide
- Pre-marked cutting guides for easy and accurate application
- Clear visibility for easy alignment over the foam

Standard Dressings

Thick and thin foam kits include green foam, port pad and transparent film



- **Thick Foam**

- Dimensions:
 - Small 7.5cm x 10cm x 3.3cm
 - Medium 12.5cm x 18cm x 3.3cm
 - Large 15cm x 26cm x 3.3cm



- **Thin Foam**

- Dimensions:
 - Small 6cm x 6cm x 2.5cm
 - Medium 15cm x 10cm x 2cm
 - Large 15cm x 23cm x 2.5cm
 - Xlarge 30cm x 60cm x 1.5cm



- **Hybrid Dressings**

- Kits include:
 - 1 Dressing
 - 1 Port pad
 - 4 Transparent Film Strips
- Hybrid Dressing Dimensions:
 - 10cm x 20cm
 - 10cm x 30cm
 - 10cm x 40cm
 - 10cm x 50cm
 - 15cm x 20cm
 - 15cm x 30cm



- **Open Abdominal Dressing**

- Kits include:
 - 1 Visceral Protection Layer
 - 2 Perforated Green Foam
 - 4 Transparent Film
 - 2 Port Pad
 - 1 Y-Connector

Specialty Dressings



- **PVA White Foam**

- Polyvinyl alcohol (PVA) foam moistened with sterile water
- Available in 3 sizes
- Dimensions:
 - Small 7.5cm x 10cm x 1cm
 - Medium 15cm x 1 cm x 1cm
 - Large 20cm x 15cm x 1cm



- **Bridge Dressing**

- Bridges dressing away from pressure areas
- Eliminates pressure at the wound site
- Dimensions:
 - 55cm length



- **Spiro Foam**

- Simple to size and allows for easier bridging
- Dimensions:
 - 15.25cm x 12.7cm x 1.5cm



- **Silver Contact Layer**


- Antimicrobial properties to prevent infection in the wound
- Dimensions:
 - Small 2cm x 1.5cm
 - Medium 5cm x 4cm
 - Large 20cm x 15cm

Canisters



- **Canister**

- 5 canister sizes available to manage all exudate levels



Accessories



- **Y-Connector**

- Allows multiple wound sites to be treated with one pump



- **Carry Bag**

- Fits all canister sizes
- Clear front screen for easy visibility
- Washable and reusable
- Adjustable shoulder strap; also for use bedside



- **Bed Hanger & IV Pole Holder**



Compact Negative
Pressure Wound Therapy
System that speeds
healing time with a quiet
and discreet design.

***Getting patients back to
enjoying daily
activities.***

Genadyne Biotechnologies

16 Midland Ave., Hicksville, NY 11801 |
Tel: +1.516.487.8787 | +1.800.208.2025

www.genadyne.com

If you have questions about our products, need support or are interested in becoming a distributor, please email or call
info@genadyne.com